

Product Datasheet

Human Angiopoietin-2 ELISA Kit (DIY Antibody Pairs) (orb1098179)

Description Human Angiopoietin-2 ELISA Kit (DIY Antibody Pairs)

Reactivity Human

Range 156 pg/ml - 10,000 pg/ml

Tested Applications ELISA

Immunogen Expression system for standard: NS0; Immunogen sequence: Y19-F496


Storage Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)

Note For research use only

Sample Types cell culture supernatants, serum, plasma (heparin, EDTA) , saliva and human milk.

Uniprot ID **O15123**

Expiration Date Please enquire.



Human
Angiopoietin-
2 EZ-Set
ELISA Kit st...