


Product Datasheet

Human CEA ELISA Kit (DIY Antibody Pairs) (orb1098194)

Description	Human CEA ELISA Kit (DIY Antibody Pairs)
Reactivity	Human
Range	312 pg/ml - 20,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: K35-A685
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
Uniprot ID	P06731
Expiration Date	Please enquire.



Human CEA
EZ-Set
ELISA Kit
standard
curv...