


Product Datasheet

Human Laminin ELISA Kit (DIY Antibody Pairs) (orb1098195)

Description	Human Laminin ELISA Kit (DIY Antibody Pairs)
Reactivity	Human
Range	156 pg/ml - 10,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: from human fibroblasts
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates and serum.
Uniprot ID	A4D0S4, P24043, Q16363, Q13751, Q13753, P25391, Q16787, P11047, Q9Y6N6, P07942, P55268, O15230
Expiration Date	Please enquire.



Human
Laminin
EZ-Set
ELISA Kit
standard ...