


Product Datasheet

Human Diablo/SMAC ELISA Kit (DIY Antibody Pairs) (orb1098203)

Description	Human Diablo/SMAC ELISA Kit (DIY Antibody Pairs)
Reactivity	Human
Range	31.2 pg/ml - 2,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: A56-D239
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
Uniprot ID	Q9NR28
Expiration Date	Please enquire.



Human
Diablo/SMAC
EZ-Set
ELISA Kit
stand...