

Product Datasheet

Human Fibronectin ELISA Kit (DIY Antibody Pairs) (orb1098209)

Description Human Fibronectin ELISA Kit (DIY Antibody Pairs)

Reactivity Human

Range 156 pg/ml - 10,000 pg/ml

Tested Applications ELISA

Immunogen Expression system for standard: from plasma


Storage Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)

Note For research use only

Sample Types cell culture supernatants, serum and plasma (heparin, EDTA, citrate).

Uniprot ID **P02751**

Expiration Date Please enquire.



Human
Fibronectin
EZ-Set
ELISA Kit
stand...