

Product Datasheet

Human NGF/NGF beta ELISA Kit (DIY Antibody Pairs) (orb1098210)

Description Human NGF/NGF beta ELISA Kit (DIY Antibody Pairs)

Reactivity Human

Range 15.6 pg/ml - 1,000 pg/ml

Tested Applications ELISA

Immunogen Expression system for standard: NS0; Immunogen sequence: S122-A241

Storage Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)


Note For research use only

Sample Types cell culture supernatants and serum.

MW 26959 MW

Uniprot ID **P01138**

Expiration Date Please enquire.



Human
NGF/NGF
beta EZ-
Set ELISA
Kit stan...