


Product Datasheet

Human CD32/FCGR2b/c ELISA Kit (DIY Antibody Pairs) (orb1098212)

Description	Human CD32/FCGR2b/c ELISA Kit (DIY Antibody Pairs)
Reactivity	Human
Range	46.9 pg/ml - 3,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: A46-P217
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
Uniprot ID	P31994, P31995
Expiration Date	Please enquire.



Human
CD32/FCGR2b/c
EZ-Set ELISA
Kit sta...