


Product Datasheet

Human IL15RA ELISA Kit (DIY Antibody Pairs) (orb1098213)

Description	Human IL15RA ELISA Kit (DIY Antibody Pairs)
Reactivity	Human
Range	156 pg/ml - 10,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: I31-T172
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA).
Uniprot ID	Q13261
Expiration Date	Please enquire.



Human
IL15RA EZ-
Set ELISA
Kit
standard
C...