


Product Datasheet

Mouse CD23/FCER2 ELISA Kit (DIY Antibody Pairs) (orb1098214)

Description	Mouse CD23/FCER2 ELISA Kit (DIY Antibody Pairs)
Reactivity	Mouse
Range	156 pg/ml - 10,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: E50-P331
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
Uniprot ID	P20693
Expiration Date	Please enquire.



Mouse
CD23/FCER2
EZ-Set
ELISA Kit
standa...