


Product Datasheet

Human TNFRSF17/BCMA ELISA Kit (DIY Antibody Pairs) (orb1098218)

Description	Human TNFRSF17/BCMA ELISA Kit (DIY Antibody Pairs)
Reactivity	Human
Range	31.2 pg/ml - 2,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: M1-A54
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
Uniprot ID	Q02223
Expiration Date	Please enquire.



Human
TNFRSF17/BCMA
EZ-Set ELISA Kit
sta...