


Product Datasheet

Human TNFRSF14/HVEM ELISA Kit (DIY Antibody Pairs) (orb1098229)

Description	Human TNFRSF14/HVEM ELISA Kit (DIY Antibody Pairs)
Reactivity	Human
Range	156 pg/ml - 10,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: L39-V202
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, tissue homogenates, serum and plasma (heparin, EDTA).
MW	30392 MW
Uniprot ID	Q92956
Expiration Date	Please enquire.



Human
TNFRSF14/HVEM
EZ Set standard
curv...