


Product Datasheet

Human Elafin/PI3 ELISA Kit (DIY Antibody Pairs) (orb1098232)

Description	Human Elafin/PI3 ELISA Kit (DIY Antibody Pairs)
Reactivity	Human
Range	31.2 pg/ml - 2,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: A23-Q117
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, tissue homogenates, serum and plasma (heparin, EDTA).
MW	12270 MW
Uniprot ID	P19957
Expiration Date	Please enquire.



Human
Elafin/PI3
EZ Set
ELISA Kit
standa...