

Product Datasheet

Human Tissue Factor/F3 ELISA Kit (DIY Antibody Pairs) (orb1098234)

Description Human Tissue Factor/F3 ELISA Kit (DIY Antibody Pairs)

Reactivity Human

Range 15.6 pg/ml - 1,000 pg/ml

Tested Applications ELISA

Immunogen Expression system for standard: NS0;
Immunogen sequence: S33-E251

Storage Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)


Note For research use only

Sample Types cell culture supernatants, serum, plasma (heparin, EDTA, citrate) , Cell lysates, and urine.

MW 33068 MW

Uniprot ID [P13726](#)

Expiration Date Please enquire.



Human
Tissue
factor/F3
EZ Set
ELISA Kit
...