


Product Datasheet

Human BMP-7 ELISA Kit (DIY Antibody Pairs) (orb1098237)

Description	Human BMP-7 ELISA Kit (DIY Antibody Pairs)
Reactivity	Human
Range	31.2 pg/ml - 2,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: S293-H431
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	bone tissue, cell culture supernatants, serum, plasma (heparin, EDTA) and urine.
MW	49313 MW
Uniprot ID	P18075
Expiration Date	Please enquire.



Human
BMP-7 EZ
Set ELISA
Kit
standard
CU...