


## Product Datasheet

# Human LOX-1/OLR1 ELISA Kit (DIY Antibody Pairs) (orb1098238)

|                            |                                                                                                                                                                                   |
|----------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Human LOX-1/OLR1 ELISA Kit (DIY Antibody Pairs)                                                                                                                                   |
| <b>Reactivity</b>          | Human                                                                                                                                                                             |
| <b>Range</b>               | 31.2 pg/ml - 2,000 pg/ml                                                                                                                                                          |
| <b>Tested Applications</b> | ELISA                                                                                                                                                                             |
| <b>Immunogen</b>           | Expression system for standard: NS0;<br>Immunogen sequence: S61-Q273                                                                                                              |
| <b>Storage</b>             | Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.) |
| <b>Note</b>                | For research use only                                                                                                                                                             |
| <b>Sample Types</b>        | cell culture supernatants, cell lysates, tissue homogenates, serum and plasma (heparin, EDTA).                                                                                    |
| <b>MW</b>                  | 30959 MW                                                                                                                                                                          |
| <b>Uniprot ID</b>          | <b>P78380</b>                                                                                                                                                                     |
| <b>Expiration Date</b>     | Please enquire.                                                                                                                                                                   |



Human  
LOX-  
1/OLR1 EZ  
Set ELISA  
Kit  
standa...