


Product Datasheet

Human Galectin-1 ELISA Kit (DIY Antibody Pairs) (orb1098239)

Description	Human Galectin-1 ELISA Kit (DIY Antibody Pairs)
Reactivity	Human
Range	156 pg/ml - 10,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: A2-D135
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, tissue homogenates, serum and plasma (heparin, EDTA).
MW	14716 MW
Uniprot ID	P09382
Expiration Date	Please enquire.



Human
Galectin-1
EZ Set
ELISA Kit
standa...