


Product Datasheet

Human CCL18/PARC ELISA Kit (DIY Antibody Pairs) (orb1098245)

Description	Human CCL18/PARC ELISA Kit (DIY Antibody Pairs)
Reactivity	Human
Range	7.8 pg/ml - 500 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: A21-A89
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants and serum.
MW	9849 MW
Uniprot ID	P55774
Expiration Date	Please enquire.



Human
CCL18/PARC
EZ Set
ELISA Kit
standa...