



## Product Datasheet

### Mouse BMP-2 ELISA Kit (DIY Antibody Pairs) (orb1098250)

|                            |                                                                                                                                                                                   |
|----------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Mouse BMP-2 ELISA Kit (DIY Antibody Pairs)                                                                                                                                        |
| <b>Reactivity</b>          | Mouse                                                                                                                                                                             |
| <b>Range</b>               | 31.2 pg/ml - 2,000 pg/ml                                                                                                                                                          |
| <b>Tested Applications</b> | ELISA                                                                                                                                                                             |
| <b>Immunogen</b>           | Expression system for standard: CHO;<br>Immunogen sequence: Q283-R396                                                                                                             |
| <b>Storage</b>             | Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.) |
| <b>Note</b>                | For research use only                                                                                                                                                             |
| <b>Sample Types</b>        | bone tissue, cell culture supernatants and serum.                                                                                                                                 |
| <b>MW</b>                  | 44525 MW                                                                                                                                                                          |
| <b>Uniprot ID</b>          | <b>P21274</b>                                                                                                                                                                     |
| <b>Expiration Date</b>     | Please enquire.                                                                                                                                                                   |



Mouse  
BMP-2 EZ  
Set ELISA  
Kit  
standard  
cu...



Mouse  
BMP-2 EZ  
Set ELISA  
Kit  
standard  
cu...