


Product Datasheet

Human Angiogenin/ANG ELISA Kit (DIY Antibody Pairs) (orb1098251)

Description	Human Angiogenin/ANG ELISA Kit (DIY Antibody Pairs)
Reactivity	Human
Range	78 pg/ml - 5,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: Q25-P147
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, tissue homogenates, serum and plasma (heparin, EDTA).
MW	16550 MW
Uniprot ID	P03950
Expiration Date	Please enquire.



Human
Angiogenin/ANG
EZ Set ELISA Kit
st...