


Product Datasheet

Human CD58 ELISA Kit (orb1098259)

Description	Human CD58 ELISA Kit
Reactivity	Human
Range	312 pg/ml - 20,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: HEK293; Immunogen sequence: F29-R215
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA, citrate)
Uniprot ID	P19256
Sensitivity	20 pg/ml
Expiration Date	Please enquire.



Human
CD58
PicoKine
ELISA Kit
standard c...