


Product Datasheet

Human DNAM-1/CD226 ELISA Kit (orb1098267)

Description	Human DNAM-1/CD226 ELISA Kit
Reactivity	Human
Range	93.7 pg/ml - 6,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: E19-A254
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA)
Uniprot ID	Q15762
Sensitivity	5 pg/ml
Expiration Date	Please enquire.



Human
DNAM-
1/CD226
PicoKine
ELISA Kit
st...