


Product Datasheet

Human sCD2 ELISA Kit (orb1098270)

Description	Human sCD2 ELISA Kit
Reactivity	Human
Range	93.7 pg/ml - 6,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: HEK293; Immunogen sequence: K25-D209
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA)
Uniprot ID	P06729
Sensitivity	5 pg/ml
Expiration Date	Please enquire.



Human sCD2
PicoKine
ELISA Kit
standard c...