

Product Datasheet

Recombinant Human CIAO2B (orb1098891)

Description	Recombinant Human CIAO2B
Reactivity	Human
Tested Applications	ELISA, WB
Preservatives	Lyophilized from a 0.2 µm filtered solution in PBS with 5% trehalose, pH7.4
Storage	-20°C for 12 months as lyophilized; 2-8°C for 1 month under sterile conditions after reconstitution
Tag	N-terminal His-IF2DI Tag
Note	For research use only
Application notes	Centrifuge the vial at 10,000 rpm for 1 minute, reconstitute at 200 µg/ml in sterile distilled water by gentle pipetting 2-3 times, don't vortex
Formula	Lyophilized from a 0.2 µm filtered solution in PBS with 5% trehalose, pH ₇ .4
Purity	Greater than 90% by SDS-PAGE gel analyses
Source	E.Coli
MW	31.3 kDa
Uniprot ID	Q9Y3D0
Expiration Date	6 months from date of receipt.