



Product Datasheet

FABP2 antibody (orb10996)


Description	FABP2 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	KLH conjugated synthetic peptide derived from chicken LFABP (21-127/127 aa)
Target	liver FABP
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	14 kDa
Uniprot ID	P12710
Dilution Range	IHC-P=1:100-500, WB=1:500-2000, ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.



IHC-P
analysis of
rat kidney
using liver...



IHC-P
analysis of
human liver
carcinoma ...



IHC-P
analysis of
mouse kidney
using liv...