

## Product Datasheet

### BrdU Cell Proliferation Assay Kit (orb109970)

**Description**
**The BrdU Cell Proliferation Assay Kit**

buffers and detection re...

**Reactivity**

All

**Tested**

ELISA

**Applications**
**Storage**

Upon receipt, store kit at -20°C in a non-frost-free freezer. For long term storage, it is recommended that you aliquot and freeze the Prediluted Anti-BrdU Detector Antibody (Component 3) and 2000x Peroxidase Goat anti-Mouse IgG (Component 5) at -20 °C. Thirty (30) minutes prior to the use of each component, thaw component. Remove the Fixative/Denaturing Solution (Component 2) and place at room temperature for at least 4 hours prior to use. The Fixative/Denaturing Solution may contain slight precipitation and its color may vary between clear to light yellow. Return the Prediluted Anti-BrdU Detector Antibody (Component 3) and 2000x Peroxidase Goat anti-Mouse IgG (Component 5) to -20°C immediately after use. All other components may be stored at 4-8°C immediately after use. Special care should be taken to keep the Prediluted Anti-BrdU Detector Antibody (Component 3) and 2000x Peroxidase Goat anti-Mouse IgG (Component 5) cold by pulling out the number of aliquots needed for the test, keeping them on ice, and leaving the remaining aliquots at -20°C.

**Note**

For research use only

**Hazard**
**Information**

This product is intended FOR RESEARCH USE ONLY, and FOR TESTS IN VITRO, not for use in diagnostic or therapeutic procedures involving humans or animals.

**Expiration Date**

12 months from date of receipt.