

Product Datasheet

Chemiluminescent BrdU Cell Proliferation Assay Kit (orb109974)

Description

Chemiluminescent BrdU Cell

Reactivity

All

Tested

ELISA

Applications
Storage

Upon receipt, store kit at -20°C in a non-frost-free freezer. For long term storage, it is recommended that you aliquot and freeze the Prediluted Anti-BrdU Detector Antibody (Component 3) and 2000x Peroxidase Goat anti-Mouse IgG (Component 5) at -20°C. Thirty (30) minutes prior to the use of each component, thaw component. Remove the Fixative/Denaturing Solution (Component 2) and place at room temperature for at least 4 hours prior to use. The Fixative/Denaturing Solution may contain slight precipitation and its color may vary between clear to light yellow. Return the Prediluted Anti-BrdU Detector Antibody (Component 3) and 2000x Peroxidase Goat anti-Mouse IgG (Component 5) to -20°C immediately after use. All other components may be stored at 4-8°C immediately after use. Special care should be taken to keep the Prediluted AntiBrdU Detector Antibody (Component 3) and 2000x Peroxidase Goat anti-Mouse IgG (Component 5) cold by pulling out the number of aliquots needed for the test, keeping them on ice, and leaving the remaining aliquots at -20°C.

Note

For research use only

Application notes

Optimal concentration should be evaluated by serial dilutions.

Hazard Information

This product is intended FOR RESEARCH USE ONLY, and FOR TESTS IN VITRO, not for use in diagnostic or therapeutic procedures involving humans or animals.

Expiration Date

12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road, Cambridge, CB5 8LA, United Kingdom

Email: info@biorbyt.com | Phone: +44 (0) 1223 859-353 | Fax: +44 (0)1223 280 240

Biorbyt LLC.

68 TW Alexander Drive
Research Triangle Park
Durham, North Carolina
27709, United States

Email: info@biorbyt.com | Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558