


Product Datasheet

Laminin beta 1 antibody (orb10998)

Description	Laminin beta 1 antibody
Species/Host	Rabbit
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-Fr, IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human Laminin beta 1 (201-300/1786 aa)
Target	Laminin Beta 1
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	196 kDa
Uniprot ID	P25391
Expiration Date	12 months from date of receipt.



Western blot analysis of mouse heart lys...