

Product Datasheet

anti- CD3 gamma antibody (orb1100475)

Description	anti- CD3 gamma antibody
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	CD3g molecule, gamma(CD3-TCR complex)
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 12 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	20-25 kDa
Uniprot ID	P09693
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.