

Product Datasheet

anti- PAR2 antibody (orb1100500)

Description	anti- PAR2 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	coagulation factor II(thrombin) receptor-like 1
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 12 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	50 kDa
Uniprot ID	P55085
Dilution Range	IHC: 1:20-1:200
Expiration Date	12 months from date of receipt.