

Product Datasheet

anti- γ -H2AX antibody (orb1100614)

| | |
|----------------------------|--|
| Description | anti- γ -H2AX antibody |
| Species/Host | Rabbit |
| Reactivity | Human, Mouse, Rat |
| Conjugation | Unconjugated |
| Tested Applications | IF, IHC, WB |
| Immunogen | A synthetic phosphorylated peptide around S139 of human H2AX |
| Form/Appearance | liquid |
| Storage | PBS with 0.02% sodium azide and 50% glycerol pH 7.3 , -20°C for 12 months (Avoid repeated freeze / thaw cycles.) |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| Purity | ≥95% as determined by SDS-PAGE |
| MW | 17 kDa |
| Uniprot ID | P16104 |
| Dilution Range | WB: 1:500-1:1000; IHC: 1:20-1:200; IF: 1:20-1:100 |
| Expiration Date | 12 months from date of receipt. |