

Product Datasheet

anti- Cathepsin B antibody (orb1100663)

Description	anti- Cathepsin B antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	cathepsin B
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 12 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	31 kDa, 25 kDa
Uniprot ID	P07858
Dilution Range	WB: 1:500-1:2000; IHC: 1:100-1:400
Expiration Date	12 months from date of receipt.