

Product Datasheet

anti- ARPC5L antibody (orb1100667)

Description	anti- ARPC5L antibody
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	actin related protein 2/3 complex, subunit 5-like
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 12 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
Uniprot ID	Q9BPX5
Dilution Range	IHC: 1:50-1:200
Expiration Date	12 months from date of receipt.