

Product Datasheet

anti- Myocilin antibody (orb1100670)

Description	anti- Myocilin antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	myocilin, trabecular meshwork inducible glucocorticoid response
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
Uniprot ID	Q99972
Dilution Range	IHC: 1:50-1:200
Expiration Date	12 months from date of receipt.