


Product Datasheet

Rat Jak2 (Tyrosine-protein kinase JAK2) ELISA Kit (orb1100754)

Description	Rat Jak2 (Tyrosine-protein kinase JAK2) ELISA Kit
Reactivity	Rat
Range	78.125-5000pg/ml
Target	Jak2
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q62689
Sensitivity	46.875pg/ml
Expiration Date	Please enquire.



Rat Jak2
(Tyrosine-
protein
kinase
JAK2) ...