


Product Datasheet

Rat Dnm1l (Dynamin-1-like protein) ELISA Kit (orb1100755)

Description	Rat Dnm1l (Dynamin-1-like protein) ELISA Kit
Reactivity	Rat
Range	78.125-5000pg/ml
Target	Dnm1l
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	O35303
Sensitivity	46.875pg/ml
Expiration Date	Please enquire.



Rat
Dnm1l
(Dynamin-
1-like
protein)
ELISA...