


Product Datasheet

Rat Cd200 (OX-2 membrane glycoprotein) ELISA Kit (orb1100757)

Description	Rat Cd200 (OX-2 membrane glycoprotein) ELISA Kit
Reactivity	Rat
Range	78.125-5000pg/ml
Target	Cd200
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	P04218
Sensitivity	46.875pg/ml
Expiration Date	Please enquire.



Rat Cd200
(OX-2
membrane
glycoprotein)
E...