


Product Datasheet

Rat NLRP3(Nod Like Receptor Pyrins-3) ELISA Kit (orb1100768)

Description	Rat NLRP3(Nod Like Receptor Pyrins-3) ELISA Kit
Reactivity	Rat
Range	0.313-20ng/ml
Target	NLRP3
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	D4A523
Sensitivity	0.188ng/ml
Expiration Date	Please enquire.



Rat
NLRP3(Nod
Like
Receptor
Pyrins-3)
EL...