


Product Datasheet

Rat DGAT2(Diacylglycerol O-acyltransferase 2)ELISA Kit (orb1100770)

Description	Rat DGAT2(Diacylglycerol O-acyltransferase 2)ELISA Kit
Reactivity	Rat
Range	31.25-2000pg/ml
Target	DGAT2
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q96PD7
Sensitivity	18.75pg/ml
Expiration Date	Please enquire.



Rat
DGAT2(Diacylglycerol
O-acyltransfera...