


Product Datasheet

Rat LTB4R(Leukotriene B4 receptor 1) ELISA Kit (orb1100772)

Description	Rat LTB4R(Leukotriene B4 receptor 1) ELISA Kit
Reactivity	Rat
Range	0.313-20ng/ml
Target	LTB4R
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q9R0Q2
Sensitivity	0.188ng/ml
Expiration Date	Please enquire.



Rat
LTB4R(Leukotriene
B4 receptor 1)
ELI...