


Product Datasheet

Rat Elovl2(Elongation of very long chain fatty acids protein 2) ELISA Kit (orb1100773)

Description	Rat Elovl2(Elongation of very long chain fatty acids protein 2) ELISA Kit
Reactivity	Rat
Range	0.156-10ng/ml
Target	Elovl2
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	D4A612
Sensitivity	0.094ng/ml
Expiration Date	Please enquire.



Rat
Elovl2(Elongation
of very long
chain...