


Product Datasheet

Rat DMP1(Dentin matrix acidic phosphoprotein 1)ELISA Kit (orb1100775)

Description	Rat DMP1(Dentin matrix acidic phosphoprotein 1)ELISA Kit
Reactivity	Rat
Range	0.156-10ng/ml
Target	DMP1
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	P98193
Sensitivity	0.094ng/ml
Expiration Date	Please enquire.



Rat
DMP1(Dentin
matrix acidic
phosphopro...