


Product Datasheet

Rat Phex (Phosphate-regulating neutral endopeptidase PHEX) ELISA Kit (orb1100776)

Description	Rat Phex (Phosphate-regulating neutral endopeptidase PHEX) ELISA Kit
Reactivity	Rat
Range	0.156-10ng/ml
Target	Phex
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	035812
Sensitivity	0.094ng/ml
Expiration Date	Please enquire.



Rat Phex
(Phosphate-regulating
neutral e...