


Product Datasheet

Rat LGR-5(leucine-rich repeat-containing g-protein coupled receptor 5) ELISA Kit (orb1100778)

Description	Rat LGR-5(leucine-rich repeat-containing g-protein coupled receptor 5) ELISA Kit
Reactivity	Rat
Range	15.625-1000pg/ml
Target	Lgr5
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	D4AC13
Sensitivity	9.375pg/ml
Expiration Date	Please enquire.



Rat LGR-5(leucine-rich repeat-containing...)