


Product Datasheet

Rat NLRP4 (Nod Like Receptor Pyrins-4) ELISA Kit (orb1100779)

Description	Rat NLRP4 (Nod Like Receptor Pyrins-4) ELISA Kit
Reactivity	Rat
Range	31.25-2000pg/ml
Target	NLRP4
Storage	2-8 °C for 6 months
Note	For research use only
Sensitivity	18.75pg/ml
Expiration Date	Please enquire.



Rat
NLRP4
(Nod Like
Receptor
Pyrins-4)
E...