


Product Datasheet

Rat NR1 (Glutamate Receptor, Ionotropic, N-Methyl-D-Aspartate 1) ELISA Kit (orb1100783)

Description	Rat NR1 (Glutamate Receptor, Ionotropic, N-Methyl-D-Aspartate 1) ELISA Kit
Reactivity	Rat
Range	0.156-10ng/ml
Target	NR1
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	P35439
Sensitivity	0.094ng/ml
Expiration Date	Please enquire.



Rat NR1
(Glutamate
Receptor,
Ionotropic,...