

Product Datasheet

Rat AMPA receptor (Glutamate receptor 1) ELISA Kit (orb1100784)

Description Rat AMPA receptor (Glutamate receptor 1) ELISA Kit

Reactivity Rat

Range 78.125-5000pg/ml

Target AMPA receptor


Storage 2-8 °C for 6 months

Note For research use only

Uniprot ID **P19490**

Sensitivity 46.875pg/ml

Expiration Date Please enquire.



Rat AMPA receptor (Glutamate receptor 1)...