


Product Datasheet

Rat SSTR2 (Somatostatin receptor type 2) ELISA Kit (orb1100785)

Description	Rat SSTR2 (Somatostatin receptor type 2) ELISA Kit
Reactivity	Rat
Range	31.25-2000pg/ml
Target	SSTR2
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	P30680
Sensitivity	18.75pg/ml
Expiration Date	Please enquire.



Rat SSTR2
(Somatostatin
receptor type
2)...