

Product Datasheet

Rat CD200 R1 (OX-2 membrane glycoprotein receptor) ELISA Kit (orb1100786)

Description Rat CD200 R1 (OX-2 membrane glycoprotein receptor) ELISA Kit

Reactivity Rat

Range 31.25-2000pg/ml

Target CD200 R1


Storage 2-8 °C for 6 months

Note For research use only

Uniprot ID [Q9ES58](https://www.uniprot.org/entry/Q9ES58)

Sensitivity 18.75pg/ml

Expiration Date Please enquire.



Rat CD200 R1
(OX-2
membrane
glycoprotein...