


Product Datasheet

Rat RF-IgM (Rheumatoid Factor IgM) ELISA Kit (orb1100787)

Description	Rat RF-IgM (Rheumatoid Factor IgM) ELISA Kit
Reactivity	Rat
Range	4.688-300IU/ml
Target	RF-IgM
Storage	2-8 °C for 6 months
Note	For research use only
Application notes	Indirect ELISA Sandwich
Sensitivity	2.813IU/ml
Expiration Date	Please enquire.



Rat RF-IgM
(Rheumatoid
Factor IgM)
ELISA...