


Product Datasheet

Rat TDO (Tryptophan-2,3-dioxygenase) ELISA Kit (orb1100789)

Description	Rat TDO (Tryptophan-2,3-dioxygenase) ELISA Kit
Reactivity	Rat
Range	0.156-10ng/ml
Target	TDO
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	P21643
Sensitivity	0.094ng/ml
Expiration Date	Please enquire.



Rat TDO
(Tryptophan-
2, 3-
dioxygenase)
EL...